

Atomic Spectroscopy Detection Limits (micrograms/liter)

Table 3. Atomic Spectroscopy Detection Limits (micrograms/liter)

Elem	Flame AA	Hg/ Hydride	GFAA	ICP-OES	ICP-MS	Elem	Flame AA	Hg/ Hydride	GFAA	ICP-OES	ICP-MS
Ag	1.5		0.005	0.6	0.002	Mo	45		0.03	0.5	0.001
Al	45		0.1	1	0.005 ^a	Na	0.3		0.005	0.5	0.0003 ^c
As	150	0.03	0.05	2	0.0006 ^b	Nb	1500			1	0.0006
Au	9		0.15	1	0.0009	Nd	1500			2	0.0004
B	1000		20	1	0.003 ^c	Ni	6		0.07	0.5	0.0004 ^c
Ba	15		0.35	0.03	0.00002 ^d	Os				6	
Be	1.5		0.008	0.09	0.003	P	75000		130	4	0.1 ^a
Bi	30	0.03	0.05	1	0.0006	Pb	15		0.05	1	0.00004 ^d
Br					0.2	Pd	30		0.09	2	0.0005
C					0.8 ^e	Pr	7500			2	0.00009
Ca	1.5		0.01	0.05	0.0002 ^d	Pt	60		2.0	1	0.002
Cd	0.8		0.002	0.1	0.00009 ^d	Rd	3		0.03	5	0.0004
Ce				1.5	0.0002	Re	750			0.5	0.0003
Cl					12	Rh	6			5	0.0002
Co	9		0.15	0.2	0.0009	Ru	100		1.0	1	0.0002
Cr	3		0.004	0.2	0.0002 ^d	S				10	28 ^j
Cs	15				0.0003	Sb	45	0.15	0.05	2	0.0009
Cu	1.5		0.014	0.4	0.0002 ^c	Sc	30			0.1	0.004
Dy	50			0.5	0.0001 ^f	Se	100	0.03	0.05	4	0.0007 ^b
Er	60			0.5	0.0001	Si	90		1.0	10	0.03 ^a
Eu	30			0.2	0.00009	Sm	3000			2	0.0002
F					372	Sn	150		0.1	2	0.0005 ^a
Fe	5		0.06	0.1	0.0003 ^d	Sr	3		0.025	0.05	0.00002 ^d
Ga	75			1.5	0.0002	Ta	1500			1	0.0005
Gd	1800			0.9	0.0008 ^g	Tb	900			2	0.00004
Ge	300			1	0.001 ^h	Te	30	0.03	0.1	2	0.0008
Hf	300			0.5	0.0008	Th				2	0.0004
Hg	300	0.009	0.6	1	0.016 ⁱ	Ti	75		0.35	0.4	0.003 ⁱ
Ho	60			0.4	0.00006	Tl	15		0.1	2	0.0002
I					0.002	Tm	15			0.6	0.00006
In	30			1	0.0007	U	1500			10	0.0001
Ir	900		3.0	1	0.001	V	60		0.1	0.5	0.0005
K	3		0.005	1	0.0002 ^d	W	1500			1	0.005
La	3000			0.4	0.0009	Y	75			0.2	0.0002
Li	0.8		0.06	0.3	0.001 ^c	Yb	8			0.1	0.0002 ^m
Lu	1000			0.1	0.00005	Zn	1.5		0.02	0.2	0.0003 ^d
Mg	0.15		0.004	0.04	0.0003 ^c	Zr	450				0.0003
Mn	1.5		0.005	0.1	0.00007 ^d						

All detection limits are given in micrograms per liter and were determined using elemental standards in dilute aqueous solution. All detection limits are based on a 98% confidence level (3 standard deviations).

All atomic absorption detection limits were determined using instrumental parameters optimized for the individual element, including the use of System 2 electrodeless discharge lamps where available. Date shown were determined on an AAnalyst™ 800.

All ICP-OES (Optima 4300/5300) detection limits were obtained under simultaneous multielement conditions with the axial view of a dual-view plasma using a cyclonic spray chamber and a concentric nebulizer.

Cold-vapor mercury detection limits were determined with a FIAS-100 or FIAS-400 flow injection system with amalgamation accessory.

The detection limit without an amalgamation accessory is 0.2 µg/L with a hollow cathode lamp, 0.05 µg/L with a System 2 electrodeless discharge lamp. (The Hg detection limited with the dedicated FIMS-100 or FIMS-400 mercury analyzers is < 0.005 µg/L without an amalgamation accessory and < 0.0002 µg/L with an amalgamation accessory.)Hydride detection limits shown were determined using anMHS-15Mercury/Hydride system.

GFAA detection limits were determined on an AAnalyst 800 using 50 µL sample volumes, an integrated and full STPF conditions.

Graphite furnace detection limits can be further enhanced by the use of replications.

Unless otherwise noted, ICP-MS detection limits were determined using an ELAN 6100/9000 equipped with Rytan™ spray chamber, Type II Cross-Flow nebulizer and nickel cones. All detection limits were determined using 3-second integration times and a minimum of 8 measurements.

Letters following an ICP-MS detection limit value refer to the use of specialize conditions or a different model instrument as follows:

a Run on ELAN DRC in standard mode using Pt cones and quartz sample introduction system. **b** Run on ELAN DRC in DRC mode using Pt cones and quartz sample introduction system. **c** Run on ELAN DRC in standard mode in Class 100 Clean Room using Pt cones and quartz sample introduction system. **d** Run Run on ELAN DRC in DRC mode in Class 100 Clean Room using Pt cones and quartz sample introduction system. **e** Using C-13. **f** Using Dy-163. **g** Using Gd-157. **h** Using Ge-74. **i** Using Hg-202. **j** Using S-34. **k** Using Te-125. **l** Using Ti-49. **m** Using Yb-173.